


# Joan's Ice Drop Spinoff

- 2 Tatting shuttles
- 1 Flat Bead (20 mm)
- Picot gauge (optional)
- ds – double stitch, p – picot, J – join
- Thread – Size 20 Lizbeth, Terra Cotta, Art. 20-673



Copyright ©2020

Inspired by [Jennifer William's Laura Heart Ice Drop](#), I created this drop that uses a slightly smaller diameter flat bead (20 mm dia. Instead of 25 mm).

For the ice drop center, work in one pass. Tat 2 heart rings then chain\* 3 times until 7 hearts. Insert the bead and tat last heart and chain to encase the bead. In the heart ring, DP is for a dropped picot – 4 of the first half of a tatted stitch, then picot followed by 4 of the second half of a tatted stitch.

\*The chains that connect the hearts are **5ds 3x (P 1ds) 5ds**.

For the row of outer chains, alternate the segments: 2 chains are with no picots (**8ds J 12ds J 8ds**) and 2 have picots (**5ds 3x (P 1ds) 5ds J 12ds J 5ds 3x (P 1ds) 5ds**). Created by Joan Thomas, [joanmariethomas@yahoo.com](mailto:joanmariethomas@yahoo.com)

For patterns with detailed instructions on how to create an ice drop, check out the [Ice Drops Addicts](#) Facebook group: <https://www.facebook.com/groups/338680599851180/>